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ab Seite 3

## Instruction manual

starting on page 22

## Mode d'emploi

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## Handleiding

vanaf pagina 58

DE

GB

FR

NL

by **GlobalMIND®**

Z 08382\_V2

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**GlobalMIND®**  
consumer electronics GmbH

11.2013 GlobalMIND GmbH, Ernst-Mantius-Str.11, 21029 Hamburg – Germany

Congratulations on deciding to purchase the VenenWalker® »basic«!

The VenenWalker® »basic« is a certified medical device which can help you prevent venous diseases and reduce the risk of spider veins and varicose veins. If you have a medical indication, please consult your GP before use!

If you have any questions regarding the VenenWalker® »basic«, contact our website: [www.venenwalker.com](http://www.venenwalker.com)

We wish you every success when using the VenenWalker® »basic«.

Regards,  
GlobalMIND!

Before using the device for the first time, please read the instructions through carefully and store them in a safe place. If the device is given to someone else, it should always be accompanied by these instructions. The manufacturer and importer do not accept any liability if the information in these instructions is not complied with!

We reserve the right to modify the device, packaging or enclosed documentation at any time in conjunction with our policy of continuous development.

## Meaning of the Symbols in These Instructions

Safety instructions are identified with this symbol. Please read them carefully and comply with the safety instructions in order to prevent personal injury and damage to property.

Tips, notes and recommendations are marked with this symbol.

Symbol for devices with protection class III. Devices with this protection class operate with safety extra-low voltage.

# Contents

Intended Use	24
Safety Instructions	24
The Proper Handling of Batteries	24
Danger of Injury	24
The VenenWalker® » basic « from GlobalMIND	25
The History of IPC (Intermittent Pneumatic Compression)	26
The IPC principle	27
Benefits of the IPC	28
Applications	29
Contraindications	30
Items Supplied and Device Overview	31
Overview of Control Unit	32
Putting Into Operation	33
Inserting/ Changing Batteries	33
Connecting the Mains Adapter	33
Use	34
After Use	35
Cleaning and Care	35
Additional Information on the VenenWalker® » basic «	36
Technical Data	38
Other Information	38
Disposal	39

## Intended Use

- ☐ The VenenWalker® »basic« must only be used for unblocking congestion in the foot and lower leg area. If you have a medical indication, please consult your GP before use!
- ☐ Use the device only in dry indoor rooms.
- ☐ Keep the device away from moist and wet areas.

## Safety Instructions

### The Proper Handling of Batteries

- ☐ Before inserting the batteries, check whether the contacts in the device and on the batteries are clean; clean them if necessary.
- ☐ Use only the battery type specified in the technical data.
- ☐ Always replace all the batteries at the same time. Do not use batteries of different types or brands, or batteries with different capacities. When replacing batteries, note the correct polarity (+/-).
- ☐ Take the batteries out of the device when they are discharged or if the device is not going to be used for a prolonged period. This will prevent any damage that might occur if a battery leaks.
- ☐ Remove dead batteries from the device immediately. Clean the contacts before you insert a new battery.
- ☐ The batteries must not be recharged or reactivated by other means, dismantled, thrown into a fire, immersed in fluids, or short-circuited.

### Danger of Injury

- ☐ Keep the device and packaging material away from children and animals. There is a risk of injury and suffocation!
- ☐ This device should not be used by people (including children) with impaired sensory or mental faculties or who do not have adequate experience and/or knowledge of such devices, unless they are supervised by a person responsible for their safety, or have received instructions from such a person on how to use the device. Children must be supervised to ensure that they do not play with the device.
- ☐ The device should not be used if it is damaged.
- ☐ Use only accessories that are supplied with the device or recommended by the original manufacturer!

## The VenenWalker® »basic« from GlobalMIND

The comprehensive massage device for the prevention of venous complaints - at a very low price.

Too much sitting, too much standing, poor nutrition: According to consistent estimates, 20 % of all Europeans have venous leg complaints and more than half is affected by advanced vein problems. This global lifestyle disease starts off harmlessly with a tingling sensation in the legs. The legs then become heavy, swollen and tired and spider veins and varicose veins appear. The entire cardiovascular system is in danger and there is a risk of ulcers, thrombosis, a pulmonary embolism and an inability to work.

How does a venous disorder occur? Veins collect the used blood in the body's tissue and transport it back to the heart. Leg muscles and venous valves tense up and relax, thus acting like a pump to transport the blood forwards into the veins. However, if a person stands or sits for a long time, or if excess strain due to surplus weight leads to congestion, the blood cannot flow away properly and the venous walls become overstretched. Swollen legs are the result. The swollen leg forms an oedema due to the fact that liquid and protein pass through the slackened venous wall and into the tissue. A doctor should be consulted as a matter of urgency at the latest when such a venous congestion becomes a continual complaint. Medicines are not able to prevent the congestion, nor reverse the abnormal slackening of the veins. Surgical stockings are a tried and tested means of supporting the weakened tissue from the exterior and preventing swollen legs. However, they should be worn at all times. "Vein walking", a special type of sport which is very popular in the USA, specifically stimulates the circulation in the veins and trains all important muscles. Again, this sport is only effective if practised consistently.

Special massages using alternating-compression devices have proven to be particularly effective and time-saving in the treatment of venous complaints. The VenenWalker® »basic« is the first ever comprehensively effective massage device which can efficiently prevent venous complaints - and this at a surprisingly low price.

The VenenWalker® »basic« comprises a control device and two double-walled leg cuffs which also cover the entire feet including the soles. The legs are lightly massaged by the VenenWalker® »basic« using an inflating and deflating air cushion which simulates the leg muscle pumping effect experienced during walking exercises. Venous reflux is improved, collected fluid is mobilised and can be flushed out.

The VenenWalker® »basic« was developed in close collaboration with medical specialists (phlebotomists), is easy to use and, if used properly, is absolutely safe and without any undesirable side effects. The entire kit contains the massage device and two leg cuffs including this detailed and easy-to-understand user manual.

## A Specialist Comments

Dr. med. Frank Sömmer, general practitioner and phlebotomist (doctor with specialist knowledge in venous disorders), with many years of experience in the treatment and prevention of venous disorders and Medical Director of the Institute for Compression Therapy in Lubeck, had the following to say about the GlobalMIND VenenWalker®:

"The VenenWalker Basic with short leg cuffs is very easy to use and can lead to a good increase in the venous return movement from the lower legs as it stimulates the important calf muscle pump."

## The History of IPC (Intermittent Pneumatic Compression)

The principle of intermittent pneumatic compression can be traced back to "quick-silver massage". In 1899, a report was already circulating about a special treatment for lymphoedemas in the upper extremities (arms). Here, the extremities were immersed in an iron cylinder filled with quicksilver.

Back in 1917, pressure increases were used in a pneumatic chamber to "expel stagnant vein blood". In the year 1929, "...intermittently inflating and deflating air cushions" were used for the first time in the treatment in bedridden patients with trophic ulcers. At the University Clinic of Leipzig, chronic venous circulation disorders with indurations (hardenings) and ulcers have been treated with the rhythmic alteration of the arterial blood flow since the beginning of the 1930s.

At around the same time, attempts were also being made in the USA to treat the consequences of artery occlusive diseases using intermittent compression.

In parallel with an increasing clinical use of the IPC, equivalent devices were further developed on a technical level right up to today's multi-chamber system for the treatment of severe lymphoedemas.

## The IPC principle

Intermittent (alternating) compression, which is not generated through manual massage but rather by a device, is known in medicine as intermittent pneumatic compression, or IPC. It uses the principle of pneumatic alternating pressures. It is prescribed by doctors for the treatment of venous or lymphatic diseases. When used to prevent thrombosis and for sports injuries, it does not carry the risks found in medication-based treatment.

The IPC simulates the function of the natural muscle pumps found in the human body. Air is pumped with alternating pressure into double-walled leg cuffs for example. The pressure is increased and decreased in defined time segments. It should be possible to adjust the pressure level to suit the user and the pressure must never exceed 100 mmHg. The treatment duration and frequency depend on the person's condition when used for cosmetic purposes, or as recommended by the doctor when used for medical reasons. However, the treatment should last for at least 20 minutes. Any medical use should, of course, be supervised by a doctor.

The shape and type of cuff used is important. Due to the proven importance of the draining of the network of veins in the sole, the cuff should cover both the entire leg and the foot.

The rhythmic increase and decrease in pressure on the tissue leads to water and protein being pumped out of the interstitial tissue. If the patient has a medical indication, a compression dressing should be applied or a surgical stocking should be worn after using the IPC as the residual protein causes blood serum to flow from the capillaries into the tissue after an IPC treatment session. This measure prevents this flow occurring. These measures do not need to be taken if you simply remain lying down after treatment. When using the equipment at home it therefore makes sense to use it before bedtime.

## Benefits of the IPC

- Improvement of venous and lymphatic macrocirculation through the passive activation of the natural muscle pump
- Improvement of the microcirculation of the skin, connective tissue and muscles
- Improvement of blood flow characteristics
- Positive flow-mechanical effects
- Increase in the blood anticoagulative protein C
- Considerable pain improvement
- Considerable boost to the wound healing process
- Quicker analgesia in cases regarding tissue fluid retention (oedemas) following operations - you can be on your feet quicker and discharged earlier
- Can promote the healing of lower-leg ulcers
- Favourable influence on lymphatic fistulas
- Carries a lower risk than other methods
- Pain-free, comfortable, uncomplicated and excellent tolerability

# Applications

## For Cosmetic Purposes

- Leg swelling caused by lack of movement, for example due standing for long periods of time in the case of sales personnel, due to sitting for longer periods for office work, or due to long-haul travel.
- The unnatural reproduction of connective tissue (fibrositis)
- Swollen legs through bedriddenness or \*pregnancy
- For varicose veins
- For muscle ache

## For medical purposes (prophylaxis)

- Swollen legs following a medical diagnosis
- Prevention of venous complaints such as varicose veins
- Prevention of varicose veins due to pregnancy
- Prevention of embolisms due to bedriddenness
- Lymphatic system swelling (lymphoedema)
- Congestion due to joint injury or overexertion, e.g. haematoma through sport injury
- Swelling following a pulled muscle or sprain
- Lower-leg ulcers caused by a venous circulation disorder
- To assist the wound healing process after an operation

- \* When using this equipment during pregnancy please note the following:  
If you notice an increased pulse frequency or experience difficulty breathing, you must stop using the device. If this happens, please contact your doctor.

## Contraindications

Do not use the IPC - or only use it only after consulting your doctor - if you have:

- Acute inflammations
- Open injuries
- Cardiac output problems with fluid retention in your legs
- Advanced arterial disease
- Suspected thrombosis
- Malignant tumors
- Pulmonary oedemas
- High blood pressure

## Items Supplied and Device Overview

- 1 2x leg cuffs
- 2 Control unit
- 3 Storage bag
- 4 2x air tubes
- 5 Mains adapter

Not shown:

- 4x 1.5 volt demo batteries, AA type

## Overview of Control Unit

- 1 Air tube connection for inflating the leg cuffs
- 2 Pressure adjuster
- 3 On/Off switch
- 4 On/Off control lamp
- 5 Air tube connection for deflating the leg cuffs

## Putting Into Operation

The device can be operated either with batteries or with the mains adapter supplied. Before use, insert the batteries or connect the mains adapter.

### Inserting/Changing Batteries

The battery compartment is on the rear of the control unit. You will need four 1.5 volt AA batteries. Use only standard batteries, never rechargeable batteries!

Four batteries are already included in this device at the time of delivery. The demo batteries supplied only have enough power to operate the device for approx. 20 minutes. Remove the transport securing device before first use.

1. Open the battery compartment.
2. Insert four 1.5 volt AA type batteries. Ensure the correct polarity (+/-) when inserting the batteries.
3. Then close the battery compartment.

Your VenenWalker® »basic« is now ready for use!

### Connecting the Mains Adapter

1. Connect the mains adapter to the control unit. The connection is located below the area for opening the battery compartment.
2. Connect the mains adapter to a properly installed socket which corresponds to the technical data of the device (see the "Technical Data" chapter).

Your VenenWalker® »basic« is now ready for use!

## Use

### CAUTION!

- ☐ To prevent damaging the air chambers, do not bend the leg cuffs! Carefully fold them after use before storing.
- ☐ Do not open the zips on the leg cuffs during use as otherwise the cuffs could get damaged.
- ☐ Please make sure that the device is always free and is not covered by a blanket or other objects.

The device can be operated both with batteries and with the supplied mains adapter. To connect the mains adapter or insert the batteries, read the chapter "Putting into Operation".

1. Put on the leg cuffs with the connections for the air tubes pointing outwards.
2. Connect the air tubes to the connections on the leg cuffs. Make sure that the tubes are connected properly and cannot slip out.
3. Connect the plug for the air tubes to the connection on the control unit for inflating the leg cuffs.
4. Set the pressure adjuster to its maximum setting.
5. Switch on the device by pressing the On/Off button. The On/Off control lamp lights up as soon as the device is turned on.
6. The leg cuffs are now alternately pumped up. If the pressure on your leg is too strong or if it feels unpleasant, set the pressure adjuster to a lower setting. The ideal pressure setting should feel strong but not unpleasant.

## How it works:

Your legs are massaged by the rhythmic and alternating inflating and deflating air cushions. This massaging of your legs has the same effect on your veins and muscles as a slightly strenuous walk. For the best results, your legs should be positioned horizontally or raised during use so that the blood can circulate unimpeded.

We recommend using the device for approximately 20 - 30 minutes.

Please note that your legs are massaged one after the other by the inflating and deflating air cushions, not simultaneously. The pressure is built up firstly in one leg cuff and then in the other. This is not a sign that the device is faulty!

## After Use

Once you have finished use, the air must be removed from the leg cuffs.

1. Set the pressure control to the highest level.
2. Switch off the device by pressing the On/Off button. The control lamp goes out.
3. Unplug the air-tube connection plug from the top socket on the control unit.
4. Now connect the plug for the air tubes to the connection on the side of the control unit.
5. Switch on the device by pressing the On/Off button. The On/Off control lamp lights up as soon as the device is turned on. The air is now pumped out of the leg cuffs.
6. Once the air has been pumped out, open the zip and remove the leg cuffs.
7. Switch off the device by pressing the On/Off button.
8. If you used the mains adapter, unplug it from both the mains socket and the device.

## Cleaning and Care

### CAUTION!

- ☐ Keep the device in a place which is protected from direct sunlight and extreme temperature fluctuations.
- ☐ Do not use any caustic or abrasive cleaning agents to clean the device. These could damage the device.

Before cleaning, always remove the air tubes from the leg cuffs and from the control unit. The leg cuffs can be cleaned with a slightly damp cloth and some household

detergent. If the leg cuffs are used for more than one person and if they come into direct contact with the skin, you should disinfect the cuffs with a standard disinfectant to avoid cross contamination. It is imperative that you make sure that no moisture enters the air inlet and outlet openings.

When used for a long time, the ends of the tubes may expand and air could leak out. In such case, you can cut off a short piece of the ends.

Only lightly fold the cuffs together. Bending the cuffs with too much force can lead to cracks.

We recommend that you only clean the control unit when it is absolutely necessary! Clean the control unit with a slightly damp cloth and, if necessary, a little household detergent.

## **Additional Information on the VenenWalker® »basic«**

### **Benefits**

- Problem-free and low-risk option of supporting decongestion treatment with compression stocking dressings
- Easy to use
- Cuffs stretching from the foot to the knee are rhythmically and alternately inflated and deflated
- The air is automatically removed from the cuffs when the plugs are switched over
- The leg is gently massaged
- The function of the leg's muscle pumps is simulated
- Not only for cosmetic problems or sports injuries
- Also for older or overweight patients
- Can be used comfortably at home on a daily basis
- Can also be used on holiday, after sitting for a long time
- Continuous acceleration of the venous outward flow for several hours
- Swellings, pains and the feeling of heaviness in the legs are reduced

Use at home for the prevention of

**CAUTION!**

☐ If you have any serious illnesses, please consult your GP before use!

- Sports injuries
- Muscle ache
- Contusions
- Sprains
- Haematomas
- Spider veins
- Varicose veins
- Simple venous congestion
- Swelling
- Feeling of heaviness
- Slight disturbances in the arterial blood flow

The intermittent compression is well proven for:

Venous, arterial, lymphatic and connective-tissue diseases

What are the results?

- Swelling is reduced
- Congestion pains subside
- The pressure in the varicose veins/spider veins is reduced, the width is reduced
- The venous valves once again close better
- To be used as a supplement to treatment for lymphoedema

## Technical Data

Cuff pressure:	Maximum 120 mmHg
Voltage supply	
In mains power mode:	
Input:	100 – 240 volts ~
Output:	6 V DC
Mains frequency:	50 Hz
Current consumption:	0.1 A
In battery mode:	
	6 volt DC (4x 1.5 volt batteries, AA type)
Temperature range	
– Operation:	15 – 40 °C
– Storage and transportation:	0 – 40 °C
Air humidity	
– Operation:	30 – 75 %
– Storage and transportation:	10 – 80 %
Dimensions:	80 mm (l) x 45 mm (w) x 190 mm (h)
Weight:	Approx. 460 grammes
Degree of protection:	Type B 
Protection class:	III 

**CE** 0123

With the CE mark, the manufacturer confirms the conformity of the product with the applicable EC directive and the compliance with the "general requirements" stipulated therein. The number 0123 alludes to the named position (certifier). Both together certify that the German Medical Devices Act (MPG) has been observed by the distributor.

## Other Information

If you have any questions regarding the VenenWalker® »basic«, contact our website: [www.venenwalker.com](http://www.venenwalker.com)

## Disposal

The packaging material is recyclable. Please dispose of it in an environmentally friendly manner and take it to a recycling centre.

Dispose of the device in an environmentally friendly manner. It should not be disposed of in the normal domestic waste. Dispose of it at a recycling centre for used electrical and electronic devices. You can obtain more information from your local authorities.

Non-rechargeable and rechargeable batteries should be removed prior to disposal of the device and disposed of separately from the device. To protect the environment, non-rechargeable and rechargeable batteries may not be disposed of with normal household waste, but must be taken to suitable collection points. Also note the relevant statutory regulations governing the disposal of batteries.

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11.2013 GlobalMIND GmbH, Ernst-Mantius-Str. 11, 21029 Hamburg – Germany